



HT4000L THE HPLC AUTOSAMPLER LAB-AUTOMATION READY

HTA - SAMPLING FOR SCIENCE | www.hta-it.com



1-800-267-8103 • instrum@chromspec.com



KEY FEATURES:

- Easy To Use
- Large sample capacity
- Application versatility
- Modular and scalable
- Lab automation ready









EASY TO USE

Easy to use, easy to install and easy to service. We believe in a "robust by design" approach: **HT4000L** – based on cutting edge technologies and DFX design – well incarnates this philosophy making itself the ideal instrument for installation in any scenario.

Full loop is performed with very high accuracy values, while **partial loop** also allows sample saving. The variable fill speed and the viscosity delay allow a wide range of sample viscosities; programmable sampling depth assures high sample recovery capacity.

analysis start, sample loading and extra wash can be managed by its dedicated keypad.

Set your method once using the complimentary HTA Autosampler Manager, then you just need to load samples and press START: experience a seamless integration with your Chromatography Data System.

LARGE SAMPLE CAPACITY

HT4000L can handle up to 176 2mlvials per run or it can work with different kinds of sample vessels. **Inject** samples out of vials (2/4/6/10/20/40 ml), test tubes, tube blocks or plates.

Save on consumables avoiding needs for sample reformatting. Choose from standard racks or enquire for a custom design: custom racks pay for themselves All the routine operations such as in a few weeks. Furthermore, racks APPLICATION VERSATILITY are easily interchangeable in order to allow quick switching across multiple configurations.

No compromise in terms of solvent capacity and solvent selection: we have a wide capacity to accommodate the most challenging requests. Solvents can be located in vials (for accurate micro- or nano-volume dispensing) or in bottles (for large-volume dispensing). Additional solvents modules may also be implemented in the field in case requirements change over the time.

No need for solvent degassing. While some autosamplers require a devoted degasser line or to degas offline wash solvent on each bottle refill, we do not. Enjoy HT4000L: less system requirement, less pain, less troubleshooting!

HT4000L fits all liquid chromatography systems. It is able to look after the widest range of applications including analytical, prep, UHPLC, IC, nano, LC/ MS... and many more!

the autosampler is equipped with an automatic 6-port 2-way injection valve for analytical applications. However, a wide selection of valves for inert or extended pressure applications are also available as per special injection **ports** for applications where no valve is needed.

MODULAR AND SCALABLE

A wide range of accessories can be sourced and added in the field at any time; they are available to implement your autosampler with. Contact your sales representative about your application requirement.

Too many options? Don't worry; you don't pay for features you don't need. We provide your LC autosampler at the best possible price; only at the time Within the standard configuration your needs change you will pay for the optional accessories to perform the new tasks without the hassle of buying a new autosampler. Your investment is future proof!

LAB AUTOMATION READY

HT4000L can fully automate your sample prep. It can work as a standard HPLC autosampler. However you may want to delegate some sample preparation tasks: HT4000L can take care of them as no standard autosampler has been able to assist.

Dilution, precolumn derivatisation, pH spiking and internal standard addition are available within the standard

UHPLC READY

UHPLC applications require an extended pressure range and a reduced amount of injected sample. **HT4000L is** optimised for such requirements and can be supplied with an UHPLC injection valve or can be upgraded in the field at a later time.

Also in this configuration **HT4000L** features extraordinary analytical performance with the highest reliability.

LARGE-VOLUME CAPABILITY

The only real autosampler choice in the market for largevolume applications.

We provide you with everything you need for your largevolume applications: large sample vessel support, large injections up to 10ml, large solvent capacity, large bore size injection valve (or injection port) for real and complete a very short time: you don't need to make a decision about automation.

BIOCOMPATIBLE

To handle more complex samples and to assure covering a wider pH range, inert flow path valves are available. No need to compromise when it comes to running metal-sensitive samples or corrosive solvents/buffers.

Furthermore, biocompatible applications can take advantage of the **double wash step capability**: as pre- and post-wash solvent in addition to A, B... F solvents you can also choose a combination of A+B, A+C.... F+E solvents for superior analytical performance (carryover adverse!).

COOLING

The Cooling Module offers the capability to maintain samples at a controlled temperature.

The cooling option upgrade can be performed in the field in cooling at the time the autosampler is sourced!





version: no extra cost required. These operations can be performed in place or in an empty vial by sample dispatching.

If standard sample preparation capabilities are not enough, get the HTAPREP option. Using HTAPREP lab automation software - you can add support to complex and custom autosampler methods. Design your own methods by Drag&Drop of autosampler functions (e.g. draw, dispense, mix, inject...) and start processing samples with the smart processing engine that allows you to optimise productivity and reagent consumption.

Endless possibilities: any sample preparation procedure or complex sample introduction become available: check the brochure of HT4000A (off-line version of HT4000L) to discover more.

TECHNICAL SPECIFICATIONS

General features

Svringe volume

Tray capacity (standard):

Tray capacity (cooled):

Maintenance: Electrical control:

Analitycal performance³ RSD

Linearity (r): Carry over2:

Nano-Volume Opt.: 5, 10, 25, 100µl 176 vials, 2ml (1 removable rack) 96 vials, 4ml (1 removable rack) 54 vials, 6/10/20ml (1 removable rack) Custom trays available 121 vials, 2ml (1 removable rack) 42 vials, 6/10/20ml (1 removable rack) Preventive counters available I AN and TTL

Standard: 100, 250, 500ul and 1, 2, 5, 10ml

Full loop: ≤ 0.25% (typically 0.15%) Partial loop: ≤ 0.5% ≥ 0.9999 $\leq 0.05\%$ (with standard wash)

Wash (Pre & Post Injection)

Washing cycle: Reagents: Stator wash:

Filling

Sample volume: Air volume: Filling speed: Viscosity delay: Bubble elimination:

Injection

Injection mode: Flush Volume Injection speed: Analysis time: Wait time between consecutive injections: typically ≤ 0.01% (with int/ext needle wash)

up to 15 (Mode: every injection/sample/step) up to 2 different reagents for wash function programmable

as low as step of 1µl³ as low as step of $1\mu l^3$ 1-100µl /sec3 0-255 up to 15 pull up strokes

Full Loop, Partial Loop

programmable

1-500µl /sec3

as low as 12sec⁴

0-999min

Diluition, derivatization, pH spiking, IS addiction Sample dispatch: step of $1\mu l^3$ step of $1\mu l^3$; up to 3 different reagents

programmable

310x705x690mm

from 4°C to 40°C

25°C +/-10°C about 5kg

0-100 minutes (0.1 minutes step)

100-240±10%Vac; 50-60Hz; 55W

+/-1°C (typically +/-0.5°C) between the vials

about 20°C below the RT (on the vial)

100-240±10%Vac; 50-60Hz; 120W

up to 15

15kg

Peltier

ves

Reagent addiction: Mixing cycle: Intermediate wash cycle: Reaction time:

Physical features

Dimensions (WxHxD): Weight: Power supply:

Cooling Module

Technology Temperature range conditioning: Temperature accuracy: Cooling capacity: Working ambient temperature: Weight Power supply: Field installation:

Additional modules:

Among the most popular modules:

- Vortex
- Barcode reader (BCR)
- Vial Heater/ Reaction Cell Gripper for vial displacement
- Internal/external needle washing
- Free software "HTA Autosampler Manager (Standard Version)" is provided for PC control; upgrade to "HTA Autosampler Manager (CFR 21 Part 11 Version)" should be purchased separately. PC is required for setup, service, method and sequence editing.

¹For volume $\geq 5\mu$ l; the values provided are for 100-1000 μ l syringes ²Tested with Uracile 100ppm in mobile phase 60% methanol : 40% water; Column C18; detector: UV fixed length 3 It depends on syringe volume; the value povided is for 100 μ syringe 4 10 μ l Partial loop injection, with Prep Ahead Mode enabled

Leading automation provider for the scientific instruments industry. HTA supplies a wide range of analyzer front-ends and sample preparation automated devices for analytical chemistry (chromatography), life sciences and clinical applications. Among its most popular products are the preparative workstations, GC and HPLC autosamplers that are commercialized worldwide thought its reseller network. HTA's quality management system is certified UNI EN ISO 9001:2008.



CHROMATOGRAPHIC

1-800-267-8103 • instrum@chromspec.com