


Disposable Sampler (UTM)



- ① 5mL Normal Cap
- ② 5mL Capture Cap
- ③ 10mL Normal Cap
- ④ Oropharyngeal Swab
- ⑤ Capture Cap Swab
- ⑥ Nasopharyngeal Swabs

Features

- Simplify specimen collection, transport and processing;
- Used for collection, storage and transportation of viruses, bacteria, chlamydiae, mycoplasma and ureaplasma. The specimen transported in the universal transport media can be used in the laboratory to perform pathogen isolation and detection;
- Disposable use.

Storage and Transport Temperature Requirement

This product could be transported at ambient temperature. Ensure samples collected are store between 2 degrees and 8 degrees Celsius, and transport to the laboratory within 48 hours. Please store the samples at minus 70 degrees Celsius if cannot transport to the laboratory within 48 hours. Please examine, inoculate or separate the sample as soon as possible after being sent to the laboratory. Those that can be inoculated or separated within 48 hours can be stored at 4 degrees Celsius. Those that cannot be inoculated or separated within 48 hours should be stored at minus 70 degrees Celsius.


Available in Canada from...

MJS
BioLynx
INC.

1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca




Quick Guide




▶▶ Collect

Collect the sample with the swab. Avoid touching the swab applicator below the breakpoint.



Oropharyngeal Swab & Capture Cap Swab:

- Gently insert swab into back of the throat and tonsillar area.
- Rub the swab over both tonsillar pillars and the posterior oropharynx, avoid touching the tongue, teeth and gums.




Nasopharyngeal Swab:

- With the patient seated, if possible, tilt the head back 70 degrees.
- Gently insert swab into the nostril along the septum floor of the nose extending straight back until the posterior nasopharynx is reached (distance from nostrils to external opening of ear).
- Rotate the swab several times while the swab is in contact with the nasopharyngeal wall.



▶▶ Snap

- Remove cap from the vial and insert the swab to the bottom of the vial.
- Holding the vial away from face, grasp the end of swab shaft and bend at a 180 degree angle to break at the breakpoint.
- Screw the cap on tightly to prevent leakage.



▶▶ Mark

Write information which you want on the vial body.



Capture Cap



- Eliminating costs associated with stocking numerous swab types.
- Just screw the cap, the swab will insert the cap automatically without human process.




Cat.No.	Description	/Pk	/Cs
202421	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Capture Cap Swabs, 10 Vials/pk, 50 Vials + 50 Swabs/cs	10	50
202422	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Capture Cap Swabs, 1 Vial +1 Swab/pk, 100 pk/cs	1	100




Normal Cap




Cat.No.	Description	/Pk	/Cs
202021	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Oropharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs	10	50
202022	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Oropharyngeal Swab, 1 Vial +1 Swab/pk, 100 pk/cs	1	100
202025	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs	10	50
202024	Disposable Sampler, 5 mL Vial with 1.2 mL UTM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 1 Vial +1 Swab/pk, 100 pk/cs	1	100
202026	Disposable Sampler, 10 mL Vial with 3.0 mL UTM, with Individually Wrapped and Sterile Oropharyngeal Swab, 1 Vial +1 Swab/pk, 100 pk/cs	1	100
202027	Disposable Sampler, 10 mL Vial with 3.0 mL UTM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 1 Vial +1 Swab/pk, 100 pk/cs	1	100
202028	Disposable Sampler, 10 mL Vial with 3.0 mL UTM, with Individually Wrapped and Sterile Oropharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs	10	50
202029	Disposable Sampler, 10 mL Vial with 3.0 mL UTM, with Individually Wrapped and Sterile Nasopharyngeal Swab, 10 Vials/pk, 50 Vials + 50 Swabs/cs	10	50




10 mL Normal Cap



5 mL Normal Cap



Oropharyngeal Swabs



Nasopharyngeal Swabs