5-second Effort ... a Perfectly-sealed Petri Dish

Shrink Seals



Fast & Easy to Use

- Slip moist band over Petri dish edge
- Each plate is 5 seconds of effort for 100% security
- No wrapping of dish required



Safe & Secure

- Tight, smooth seal every time
- No gaps or ripping
- Pathogenic microorganisms safely contained



Preserves Media Integrity

- Blocks evaporation during extended incubation
- Extends shelf-life of sensitive or pre-made plates

Other products available at Scientific Device Laboratory



- Avoid bubbling or foggy tape preps with Fungi-Tape (#745)
- Create beautiful permanent mounts with built-in stain MycoPerm (#753)
- Shrink Seals (#261)

For more information see the catalog on our website.







inquire for FREE SAMPLE of Shrink Seals




scientificdevice.com

Available in Canada from...












MycoPerm

- Integrate your permanent mycology mounting & stain into one step
- Protect preps with plasticizers for a durable, long-life with this novel stain
- Eliminate messy nail polish
 - 745-1 Mycology Mounting Kit BLUE 1 Fungitape + 12ml MycoPerm BLUE/pkg
 - 745-9 Mycology Mounting Kit RED 1 Fungitape + 12ml MycoPerm RED/pkg
 - 752 MycoPerm RED 12ml bottle/pkg
 - 753 MycoPerm BLUE 12ml bottle/pkg

The Hines Method

... an innovative technique

The Mycology Laboratory at Hines VA Medical Center in suburban Chicago just loves their MycoPerm. Their innovative technique involves one drop of MycoPerm RED with acid fuchsin and one drop of MycoPerm BLUE with Lacto Phenol Cotton blue on each slide. As they scan the permanent mount they have a spectrum of color where fungal structures are sharper and even easier to visualize.



Fungi-Tape

- Improve tape-preps for fungi slide mounts with specially formulated tape that does not disintegrate, fog or bubble
- Select colonies with easy to dispense, pliable adhesive tape just form loop
- Professionalize process with strong tape, each roll contains approximately 360 coverslips

745 Fungi-Tape - 1 roll/pkg



scientificdevice.com

