



Phthalate Standards Reference Guide



AccuStandard®

Available in Canada from ...



**CHROMATOGRAPHIC
SPECIALTIES INC.**

1-800-267-8103 • www.chromspec.com • tech@chromspec.com

Phthalate Background

Benzene dicarboxylic acid is equivalent to phthalic acid. Reacting phthalic acid with a variety of alcohols results in the synthesis of a group of chemicals designated as phthalic acid esters or phthalates.

Phthalates are used primarily as plasticizers. Plasticizers lower the glass transition temperature of a plastic/polymer and impart flexibility, durability and longevity to these types of products by acting as softening agents.

Due to their low-cost, versatility and effectiveness, phthalates are widely used in plastics manufacturing, in pharmaceutical coatings, in all types of packaging, in inks, in textiles and as gelling agents. They are end-use components of electronics, paints, adhesives, building materials, cleaning products and toys to name just a few.

Phthalates are now separated into two distinct classes according to the length of the precursor alcohol. The lower molecular weight (LMW) phthalates, including di-butyl, benzyl butyl and diethyl hexyl are made from alcohols with three to six carbon backbones.

Unfortunately, the LMW phthalates are easily released into the environment because there is no chemical bond between the phthalates and the plastic/polymer matrix. Leaching and atmospheric release of these compounds increases as the substrate ages and/or weathers; and has resulted in phthalates becoming a major environmental contaminant. This is important because phthalates are considered to be potential endocrine-disrupting agents (1). Human exposure to phthalates may be through direct contact, ingestion or inhalation. Concern over the adverse health effects has prompted regulatory changes and lead to a permanent ban of these plasticizers in baby-care products and toys (2).

Such a large-scale health concern has led to the development of analytical methods for phthalates in a variety of matrices. The majority of these methods focus on the analysis of the LMW phthalates ranging from mono/diethyl to mono/dioctyl and, in particular, dibutyl and diethyl hexyl phthalates (3). Diethyl hexyl phthalate has been the dominant plasticizer and is the largest volume phthalate in the global market. It is used as a standard for comparison for the performance of other types of plasticizers.

All of the above-mentioned phthalates are single isomer compounds which can be analyzed via straight-forward GC/MS methods yielding a single chromatographic peak for each compound.



Phthalate Contents	
Phthalates	1-3
Isophthalates	3
Terephthalates	3
Monophthalates	3-4
Deuterated Phthalates	5
Phthalates - Tech Mixes	6
EPA Methods	7
Set, Kit and Mixtures	Back

1. S. Jobling et al., *Environ. Health Perspect.*, **103** (6), 582-587 (1995)
2. *Chemical & Engineering News*, vol. 89, no 22, page 28 (May 30 2011)
3. H. Fromme et al., *Water Research*, **36** (6), 1429-1438 (2002)

The perfect companions for your analysis!

Handbook for the Chemical Analysis of Plastic and Polymer Additives, 2nd Edition



This reference handbook contains the compounds with important reference data to aid in testing and compliance. There is also information to help with real world examples, tips for analysis in challenging matrices, and much more.

The new book will be available Spring 2015



Cat. No: PLAS-CRC-BOOK2

Each Compound has:

Chemical Information

- Structure
- CAS Number (where applicable)
- RTECS Number (where available)
- Formula
- Molecular Weight
- IUPAC Name, common names and popular brand names (where available)

Physical Properties

- Appearance
- Melting and Boiling Points
- Stability
- Solubilities in several common solvents


Other Important Information

- Application
- Regulatory
- Environmental Impact
- Point of Release
- Toxicological Data

Analytical Data

- Mass Spectrum with key ions tabulated
- Chromatogram with conditions

Plastic Additive Standards Guide Includes phthalates and other plasticizers



Both the Handbook and Guide are organized into classes by additive type. Manufacturers can easily find Standards that match their particular application and product formulation for the following product categories:


- Medical Devices
- Food Packaging
- Pharmaceutical Packaging
- Toys
- Wire and Cable
- etc.

Product Group including:

Accelerants, Antidegradants, Antifoams, Antioxidants, Antiozonants, Blowing Agents, Coupling Agents, Crosslinking Agents, Flame Retardants, Plasticizers, Processing Aids, Retarders, Stearates, UV Stabilizers, Vegetable Oils, Deuterated Phthalates

Phthalates


Benzyl butyl phthalate



CAS 85-68-7 MF C₁₉H₂₀O₄ MW 312.26 PS L
SG 1.13 g/cm³ MP -35 °C BP 370-380 °C
FP 198 °C

Matrix	Cat. No.	Unit
NEAT	ALR-082N	100 mg
100 µg/mL in MeOH	ALR-082S	1 mL
5 mg/mL in MeOH	AS-E0065	1 mL


bis(2-n-Butoxyethyl)phthalate



CAS 117-83-9 MF C₂₀H₃₀O₆ MW 366.45 PS L
SG 1.06 g/cm³ MP N/A BP 270 °C FP 205 °C

Matrix	Cat. No.	Unit
NEAT	J-112	100 mg

bis(2-Ethoxyethyl)phthalate




CAS 605-54-9 MF C₁₆H₂₂O₆ MW 310.34 PS S
SG 1.12 g/cm³ MP 34 °C BP 345 °C FP N/A

Matrix	Cat. No.	Unit
NEAT	J-111	100 mg

bis(2-Ethylhexyl)phthalate (DEHP)


Di(2-Ethylhexyl) phthalate



CAS 117-81-7 MF C₂₄H₃₈O₄ MW 390.56 PS L
SG 0.98 g/cm³ MP -50 °C BP 361 °C FP 204 °C

Matrix	Cat. No.	Unit
NEAT	ALR-097N	100 mg
100 µg/mL in MeOH	ALR-097S	1 mL
1 mg/mL in MeOH	APP-9-029-10X	1 mL


bis(2-Methoxyethyl)phthalate



CAS 117-82-8 MF C₁₄H₁₈O₆ MW 282.29 PS L
SG 1.17 g/cm³ MP N/A BP 230 °C FP 121 °C

Matrix	Cat. No.	Unit
NEAT	J-106	100 mg


bis(4-Methyl-2-pentyl)phthalate



CAS 146-50-9 MF C₂₀H₃₀O₄ MW 334.45 PS L
SG 1.01 g/cm³ BP 370-380 °C FP 180 °C

Matrix	Cat. No.	Unit
NEAT	J-109	100 mg


2-Butoxy-2-oxoethyl butyl phthalate



CAS 85-70-1 MF C₁₈H₂₄O₆ MW 336.38 PS L
SG 1.10 g/cm³ MP N/A BP 345 °C FP 199 °C

Matrix	Cat. No.	Unit
NEAT	J-115	100 mg

Diallyl phthalate




CAS 131-17-9 MF C₁₄H₁₆O₄ MW 246.26 PS L
SG 1.11 g/cm³ MP 16 °C BP 165 °C FP 165 °C

Matrix	Cat. No.	Unit
NEAT	J-002	100 mg

Property Key

CAS	Chemical Abstract Service Number	SG	Specific Gravity (g/cm ³)
MF	Molecular Formula	MP	Melting Point (°C)
MW	Molecular Weight	BP	Boiling Point (°C)
PS	Physical State (Solid, Liquid)	FP	Flash Point (°C)


Diamyl phthalate



CAS 131-18-0 MF C₁₈H₂₆O₄ MW 306.40 PS L
SG 1.03 g/cm³ BP 342 °C FP 190 °C

Matrix	Cat. No.	Unit
NEAT	ALR-098N	100 mg
100 µg/mL in MeOH	ALR-098S	1 mL


Dibenzyl phthalate



CAS 523-31-9 MF C₂₂H₁₈O₄ MW 346.38 PS S
SG 1.25 g/cm³ MP 40-42 °C BP >400 °C
FP >150 °C

Matrix	Cat. No.	Unit
NEAT	J-104	100 mg


Dibutyl phthalate



CAS 84-74-2 MF C₁₆H₂₂O₄ MW 278.34 PS L
SG 1.05 g/cm³ MP -35 °C BP 337-340 °C
FP 177 °C

Matrix	Cat. No.	Unit
NEAT	J-003	100 mg
100 µg/mL in MeOH	APP-9-063	1 mL
1 mg/mL in MeOH	APP-9-063-10X	1 mL
5 mg/mL in MeOH	AS-E0066	1 mL

Dicyclohexyl phthalate




CAS 84-61-7 MF C₂₀H₂₆O₄ MW 330.42 PS S
SG 1.14 g/cm³ MP 61-66 °C BP 235 °C
FP 207 °C

Matrix	Cat. No.	Unit
NEAT	J-004	100 mg
100 µg/mL in MeOH	ALR-099S	1 mL
1 mg/mL in AcCN	AS-E0318	1 mL

Phthalates continued on next page

Phthalates (continued)


Didodecyl phthalate



CAS 2432-90-8 MF C₃₂H₅₄O₄ MW 502.77 PS L or S
SG 1.05 g/cm³ MP 21-23 °C BP N/A FP >200 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-018N	100 mg
100 µg/mL in Hexane	PHTH-018S	1 mL


Diethyl phthalate



CAS 84-66-2 MF C₁₂H₁₄O₄ MW 222.24 PS L
SG 1.12 g/cm³ MP -3 °C BP 172 °C FP 160 °C

Matrix	Cat. No.	Unit
NEAT	J-005	100 mg
100 µg/mL in MeOH	APP-9-081	1 mL
1 mg/mL in MeOH	APP-9-081-10X	1 mL
5 mg/mL in MeOH	AS-E0068	1 mL


Dihexyl phthalate



CAS 84-75-3 MF C₂₀H₃₀O₄ MW 334.45 PS L
SG 1.01 g/cm³ BP 185-187 °C FP 200 °C

Matrix	Cat. No.	Unit
NEAT	ALR-100N	100 mg
100 µg/mL in MeOH	ALR-100S	1 mL


Diphenyl phthalate



CAS 84-62-8 MF C₂₀H₁₄O₄ MW 318.32 PS S
SG 1.24 g/cm³ MP 74-76 °C BP 255 °C FP 256 °C

Matrix	Cat. No.	Unit
NEAT	J-013	100 mg


Diisopentyl phthalate **NEW**



CAS 605-50-5 MF C₁₈H₂₆O₄ MW 306.40 PS L
SG 1.03 g/cm³ FP 167 °C

Matrix	Cat. No.	Unit
NEAT	J-127	100 mg


Diisopropyl phthalate



CAS 605-45-8 MF C₁₄H₁₈O₄ MW 250.29 PS L

Matrix	Cat. No.	Unit
NEAT	PHTH-019N	100 mg
100 µg/mL in MeOH	PHTH-019S	1 mL


Dimethyl phthalate



CAS 131-11-3 MF C₁₀H₁₀O₄ MW 194.18 PS L
SG 1.19 g/cm³ MP 2-6 °C BP 282-284 °C
FP 295 °C

Matrix	Cat. No.	Unit
NEAT	J-010	100 mg
100 µg/mL in MeOH	APP-9-088	1 mL
1 mg/mL in MeOH	APP-9-088-10X	1 mL
5 mg/mL in MeOH	AS-E0069	1 mL
0.1 mg/mL in EtOAc	M-8032-IS	1 mL


Di-n-heptyl phthalate



CAS 3648-21-3 MF C₂₂H₃₄O₄ MW 362.50 PS L
SG 0.99 g/cm³ MP N/A BP 195 °C FP 113 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-020N	100 mg
100 µg/mL in MeOH	PHTH-020S	1 mL


Di-n-octyl phthalate



CAS 117-84-0 MF C₂₄H₃₈O₄ MW 390.56 PS XX
SG X.XX g/cm³ MP XXX °C FP XXX °C

Matrix	Cat. No.	Unit
NEAT	J-011	100 mg


Di-n-propyl phthalate



CAS 131-16-8 MF C₁₄H₁₈O₄ MW 250.29 PS L
SG 1.08 g/cm³ MP N/A BP 317-318 °C FP 109 °C

Matrix	Cat. No.	Unit
NEAT	J-100	100 mg


Diisobutyl phthalate



CAS 84-69-5 MF C₁₆H₂₂O₄ MW 278.34 PS L or S
SG 1.04 g/cm³ MP N/A BP 327 °C FP 109 °C

Matrix	Cat. No.	Unit
NEAT	J-113	100 mg

Diundecyl phthalate



CAS 3648-20-2 MF C₃₀H₅₀O₄ MW 474.72 PS L or S
SG 0.95 g/cm³ MP 15 °C BP 472 °C FP 239 °C


Matrix	Cat. No.	Unit
NEAT	PHTH-021N	100 mg
100 µg/mL in MeOH	PHTH-021S	1 mL

Property Key

CAS	Chemical Abstract Service Number	SG	Specific Gravity (g/cm ³)
MF	Molecular Formula	MP	Melting Point (°C)
MW	Molecular Weight	BP	Boiling Point (°C)
PS	Physical State (Solid, Liquid)	FP	Flash Point (°C)

Isophthalates


Dimethyl isophthalate



CAS 1459-93-4 **MF** C₁₀H₁₀O₄ **MW** 194.18 **PS** S
SG 1.18 g/cm³ **MP** 64-68 °C **BP** 282-285 °C
FP 148 °C

Matrix	Cat. No.	Unit
NEAT	J-009	100 mg

Diphenyl isophthalate




CAS 744-45-6 **MF** C₂₀H₁₄O₄ **MW** 318.32 **PS** S
SG 1.24 g/cm³ **MP** 136-138 °C **FP** 256 °C

Matrix	Cat. No.	Unit
NEAT	J-012	100 mg

Terephthalates


bis(2-Ethylhexyl)terephthalate



CAS 6422-86-2 **MF** C₂₄H₃₈O₄ **MW** 390.56 **PS** L
SG 0.99 g/cm³ **MP** N/A **BP** 400 °C **FP** 212 °C

Matrix	Cat. No.	Unit
NEAT	J-121	100 mg


Diethyl terephthalate



CAS 636-09-9 **MF** C₁₂H₁₄O₄ **MW** 222.24 **PS** S
SG 1.15 g/cm³ **MP** 43-47 °C **BP** 142 °C
FP >150 °C

Matrix	Cat. No.	Unit
NEAT	J-123	100 mg

Dimethyl terephthalate




CAS 120-61-6 **MF** C₁₀H₁₀O₄ **MW** 194.18 **PS** S
SG 1.36 g/cm³ **MP** 139-141 °C **BP** 288 °C
FP 151 °C

Matrix	Cat. No.	Unit
NEAT	J-101	100 mg

Monophthalates

Mono-phthalate esters are the primary phthalate metabolites formed via hydrolysis of one ester bond. It is these compounds that are thought to be toxic agents; and are receiving interest as a possible human health issue. Studies have shown that they can produce estrogenic and immune-suppressive effects in humans.


Monobenzyl phthalate (mBzP) **NEW**



CAS 2528-16-7 **MF** C₁₅H₁₂O₄ **MW** 256.25 **PS** S
SG 1.28 g/cm³ **MP** 106 °C **BP** 441 °C **FP** 168 °C

Matrix	Cat. No.	Unit
NEAT	ALR-134N	100 mg
100 µg/mL in AcCN	ALR-134S-CN	1 mL


Monobutyl phthalate (mBP)



CAS 131-70-4 **MF** C₁₂H₁₄O₄ **MW** 222.24 **PS** S
SG 1.17 g/cm³ **MP** 73 °C **BP** 350-354 °C **FP** 138 °C

Matrix	Cat. No.	Unit
NEAT	ALR-135N	100 mg
100 µg/mL in AcCN	ALR-135S-CN	1 mL

Monocyclohexyl phthalate (mBP) **NEW**



CAS 7517-36-4 **MF** C₁₄H₁₆O₄ **MW** 248.27 **PS** S
SG 1.24 g/cm³ **MP** 89-91 °C **BP** 410 °C **FP** 154 °C

Matrix	Cat. No.	Unit
NEAT	ALR-178N	100 mg
100 µg/mL in AcCN	ALR-178S-CN	1 mL


Technical Note

AccuStandard offers eight mono-phthalates including the mono-ethylhexyl (mEHP) which is the metabolite of the plasticizer with the greatest yearly production and use on a global basis.

Monophthalates continued on next page

Monophthalates (continued)


Monoethyl phthalate (mEP)



CAS 2306-33-4 **MF** C₁₀H₁₀O₄ **MW** 194.18 **PS** S
SG 1.24 g/cm³ **MP** 101 °C **BP** 329 °C **FP** 135 °C

Matrix	Cat. No.	Unit
NEAT	ALR-137N	100 mg
100 µg/mL in AcCN	ALR-137S-CN	1 mL


Monoethylhexyl phthalate (mEHP)



CAS 4376-20-9 **MF** C₁₆H₂₂O₄ **MW** 278.34 **PS** S
SG 1.09 g/cm³ **MP** 142 °C **BP** 390-395 °C
FP 144 °C

Matrix	Cat. No.	Unit
NEAT	ALR-138N	100 mg
100 µg/mL in AcCN	ALR-138S-CN	1 mL


Mono-2-heptyl phthalate **NEW**



CAS N/A **MF** C₁₅H₂₀O₄ **MW** 264.32 **PS** S

Matrix	Cat. No.	Unit
NEAT	ALR-143N	100 mg
100 µg/mL in AcCN	ALR-143S-CN	1 mL


Monoethyl phthalate **NEW**



CAS 24539-57-9 **MF** C₁₄H₁₆O₄ **MW** 250.29 **PS** S
SG 1.12 g/cm³ **MP** 133 °C **BP** 375-380 °C
FP 142 °C

Matrix	Cat. No.	Unit
NEAT	ALR-175N	100 mg
100 µg/mL in AcCN	ALR-175S-CN	1 mL


Monoisobutyl phthalate **NEW**



CAS 30833-53-5 **MF** C₁₂H₁₄O₄ **MW** 222.24 **PS** S
SG 1.17 g/cm³ **MP** 78-80 °C **BP** 356-357 °C
FP 135 °C

Matrix	Cat. No.	Unit
NEAT	ALR-176N	100 mg
100 µg/mL in AcCN	ALR-176S-CN	1 mL


Monoisononyl phthalate **NEW**



CAS N/A **MF** C₁₇H₂₄O₄ **MW** 292.37 **PS** S

Matrix	Cat. No.	Unit
NEAT	ALR-142N	100 mg
100 µg/mL in AcCN	ALR-142S-CN	1 mL


Monoisopropyl phthalate **NEW**



CAS 35118-50-4 **MF** C₁₁H₁₂O₄ **MW** 208.21 **PS** S
SG 1.20 g/cm³ **MP** 100-104 °C **BP** 343-344 °C
FP 133 °C

Matrix	Cat. No.	Unit
NEAT	ALR-179N	100 mg
100 µg/mL in AcCN	ALR-179S-CN	1 mL


Monomethyl phthalate



CAS 4376-18-5 **MF** C₉H₈O₄ **MW** 180.16 **PS** S
SG 1.29 g/cm³ **MP** 81-84 °C **BP** 315-316 °C
FP 135 °C

Matrix	Cat. No.	Unit
NEAT	ALR-139N	100 mg
100 µg/mL in AcCN	ALR-139S-CN	1 mL


Monoethyl phthalate **NEW**



CAS 5393-19-1 **MF** C₁₆H₂₂O₄ **MW** 278.34 **PS** S
SG 1.09 g/cm³ **MP** 149 °C **BP** 400 °C **FP** 146 °C

Matrix	Cat. No.	Unit
NEAT	ALR-141N	100 mg
100 µg/mL in AcCN	ALR-141S-CN	1 mL

Mono-n-pentyl phthalate **NEW**



CAS 24539-56-8 **MF** C₁₃H₁₈O₄ **MW** 236.26 **PS** S
SG 1.48 g/cm³ **MP** 125-126 °C **BP** 376-377 °C
FP 140 °C

Matrix	Cat. No.	Unit
NEAT	ALR-177N	100 mg
100 µg/mL in AcCN	ALR-177S-CN	1 mL

Property Key


CAS	Chemical Abstract Service Number	SG	Specific Gravity (g/cm ³)
MF	Molecular Formula	MP	Melting Point (°C)
MW	Molecular Weight	BP	Boiling Point (°C)
PS	Physical State (Solid, Liquid)	FP	Flash Point (°C)



Deuterated Phthalates

AccuStandard offers eleven deuterated phthalates which can be used as internal standards. To simplify the ordering process, the native and the corresponding deuterated compound are packaged as sets at a reduced price


Dibenzylphthalate-d₄



CAS 1015854-62-2 MF C₂₀H₁₄D₄O₄ MW 350.40
PS S MP 40-42 °C BP 276-278 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-001N	5 mg
100 µg/mL in MeOH	PHTH-D4-001S	1 mL


Dicyclohexyl phthalate-3,4,5,6-d₄



CAS 358731-25-6 MF C₂₀H₂₂D₄O₄ MW 334.44 PS S
SG 1.16 g/cm³ MP 65-67 °C FP 207 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-004N	5 mg
100 µg/mL in MeOH	PHTH-D4-004S	1 mL


Diethyl phthalate-3,4,5,6-d₄



CAS 93952-12-6 MF C₁₀H₁₀D₄O₄ MW 226.26 PS L
SG 1.14 g/cm³ MP -3 °C FP 160 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-005N	5 mg
100 µg/mL in MeOH	PHTH-D4-005S	1 mL


Di-iso-butyl phthalate-3,4,5,6-d₄



CAS 358730-88-8 MF C₁₆H₁₈D₄O₄ MW 282.37 PS L
MP N/A BP 327 °C FP 109 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-003N	5 mg
100 µg/mL in MeOH	PHTH-D4-003S	1 mL


Dimethyl phthalate-3,4,5,6-d₄



CAS 93951-89-4 MF C₁₀H₈D₄O₄ MW 198.21 PS L
SG 1.20 g/cm³ MP 2 °C BP 282 °C FP 147 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-007N	5 mg
100 µg/mL in MeOH	PHTH-D4-007S	1 mL


Di-n-butyl phthalate-d₄



CAS 93952-11-5 MF D₁₆H₁₈D₄O₄ MW 282.37 PS L
SG 1.07 g/cm³ MP < 25 °C BP 336 °C FP 171 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-002N	5 mg
100 µg/mL in MeOH	PHTH-D4-002S	1 mL


Di-n-hexyl phthalate-3,4,5,6-d₄



CAS 1015854-55-3 MF C₂₀H₂₆D₄O₄ MW 338.47
PS L SG 1.01 g/cm³ MP N/A BP 185-187 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-006N	5 mg
100 µg/mL in MeOH	PHTH-D4-006S	1 mL


Di-n-octyl phthalate-3,4,5,6-d₄



CAS 93952-13-7 MF C₂₆H₃₄D₄O₄ MW 394.58 PS L
SG 0.96 g/cm³ MP -50 °C BP 390 °C FP 205 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-008N	5 mg
100 µg/mL in MeOH	PHTH-D4-008S	1 mL


Di-n-pentyl phthalate-3,4,5,6-d₄



CAS 358730-89-9 MF C₁₈H₂₂D₄O₄ MW 310.42 PS L
MP N/A BP 360 °C FP >110 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-009N	5 mg
100 µg/mL in MeOH	PHTH-D4-009S	1 mL


Di-n-propyl phthalate-3,4,5,6-d₄



CAS 358731-29-0 MF C₁₄H₁₄D₄O₄ MW 254.31 PS L
SG 1.08 g/cm³ MW N/A BP 317-318 °C FP >110 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-010N	5 mg
100 µg/mL in MeOH	PHTH-D4-010S	1 mL

bis(2-Ethylhexyl) phthalate-3,4,5,6-d₄



CAS 93951-87-2 MF C₂₄H₃₄D₄O₄ MW 394.58 PS L
SG 0.98 g/cm³ MP -50 °C BP 384 °C FP 207 °C

Matrix	Cat. No.	Unit
NEAT	PHTH-D4-011N	5 mg
100 µg/mL in MeOH	PHTH-D4-011S	1 mL

Set includes 11 Deuterated Phthalates

Neat Set	PHTH-D4N-SET	11 x 5 mg
Solution Set	PHTH-D4S-SET	11 x 1 mL

Other compounds are available.
contact our Technical Service if you
require additional deuterated or other
labeled compounds.

Phthalates - Technical Mixtures

The high molecular weight (HMW) phthalates have more than six carbons in the backbone; and are synthesized from phthalic acid and mixtures of C9 and C10 alcohols. The two major HMW products are diisononyl phthalate (DINP) and diisodecyl phthalate (DIDP).

Attention has now turned to the analysis of these compounds as they are becoming major players in the plasticizer marketplace. However, due to the synthesis process, GC separation of DINP and DIDP results in a cluster of peaks corresponding to different isomers. Consequently, different analytical approaches based on soft ionization techniques and MS detection have been documented in the literature (1). These new approaches can provide another tool to scrutinize the amounts, environmental fate and potential health effects of these HMW plasticizers.

1. David, F., Sandra, P. and Hancock, P., *Current Trends in Mass Spectrometry*, May 2011)

Compound	CAS No.	Conc.	Matrix	Cat. No.
Benzyl 2-ethylhexyl phthalate	27215-22-1	100 mg	NEAT	ALR-165N
		100 µg/mL	MeOH	ALR-165S
n-Butyl benzyl phthalate NEW	85-68-7	10 mg	NEAT	PHTH-014N
		100 µg/mL	MeOH	PHTH-014S
Butyl cyclohexyl phthalate	84-64-0	100 mg	NEAT	J-122
n-Butyl isobutyl phthalate NEW	17851-53-5	10 mg	NEAT	PHTH-013N
		100 µg/mL	MeOH	PHTH-013S
Butyl octyl phthalate	84-78-6	100 mg	NEAT	J-001
Decyl octyl phthalate NEW	119-07-3	10 mg	NEAT	PHTH-012N
		100 µg/mL	MeOH	PHTH-012S
Didecyl phthalate NEW	84-77-5	100 mg	NEAT	J-120
Diisodecyl phthalate	26761-40-0	100 mg	NEAT	ALR-101N
		100 µg/mL	MeOH	ALR-101S
Diisoheptyl phthalate	71888-89-6	50 mg	NEAT	PHTH-017N
		100 µg/mL	MeOH	PHTH-017S
Diisohexyl phthalate	68515-50-4	100 mg	NEAT	J-007
Diisononyl phthalate (C8 to C10 Isomers)	68515-48-0	100 mg	NEAT	ALR-102N
		100 µg/mL	MeOH	ALR-102S
Diisooctyl phthalate (C8 Isomers)	27554-26-3	100 mg	NEAT	ALR-103N
		100 µg/mL	MeOH	ALR-103S
Dinonyl phthalate NEW	84-76-4	100 mg	NEAT	J-105
Hexyl-2-ethylhexyl phthalate NEW	75673-16-4	100 mg	NEAT	J-016
Isobutyl benzyl phthalate NEW	72170-45-7	10 mg	NEAT	PHTH-015N
		100 µg/mL	MeOH	PHTH-015S
Isobutyl cyclohexyl phthalate	5334-09-8	100 mg	NEAT	J-014
Pentyl isopentyl phthalate NEW	776297-69-9	10 mg	NEAT	PHTH-016N
		100 µg/mL	MeOH	PHTH-016S
n-Octyl n-decyl phthalate	119-07-3	100 mg	NEAT	J-015



EPA Methods - Phthalate Standards

Method 506 Phthalate Esters by PID

Phthalate Esters

M-506		1 x 1 mL
M-506-PAK	SAVE	5 x 1 mL
1.0 mg/mL each in Isooctane		
Benzyl butyl phthalate	bis(2-Ethylhexyl)adipate	
Dimethyl phthalate	bis(2-Ethylhexyl)phthalate	
Diethyl phthalate	Di- <i>n</i> -octyl phthalate	
Di- <i>n</i> -butyl phthalate		

M-506-QC		1 x 1 mL
M-506-QC-PAK	SAVE	5 x 1 mL
At stated conc. in MeOH		

Benzyl butyl phthalate (0.25 mg/mL)	bis(2-Ethylhexyl)adipate (1.2 mg/mL)
Dimethyl phthalate (0.1 mg/mL)	bis(2-Ethylhexyl)phthalate (0.25 mg/mL)
Diethyl phthalate (0.1 mg/mL)	Di- <i>n</i> -octyl phthalate (0.65 mg/mL)
Di- <i>n</i> -butyl phthalate (0.1 mg/mL)	

Method 606 Phthalate Esters by GC/ECD

M-606		1 x 1 mL
M-606-PAK	SAVE	5 x 1 mL
0.2 mg/mL each in MeOH		

Benzyl butyl phthalate	Di- <i>n</i> -butyl phthalate
Dimethyl phthalate	Di- <i>n</i> -octyl phthalate
Diethyl phthalate	bis(2-Ethylhexyl)phthalate

Method 8060 Phthalate Esters by GC/ECD

Phthalate Esters

M-8060		1 x 1 mL
M-8060-PAK	SAVE	5 x 1 mL
2.0 mg/mL each in Isooctane		
Benzyl butyl phthalate	Di- <i>n</i> -butyl phthalate	
Diethyl phthalate	Di- <i>n</i> -octyl phthalate	
Dimethyl phthalate	bis(2-Ethylhexyl)phthalate	

M-8060-QC		1 x 1 mL
M-8060-QC-PAK	SAVE	5 x 1 mL
At stated conc. in MeOH		

Benzyl butyl phthalate (0.1 mg/mL)	Di- <i>n</i> -butyl phthalate (0.25 mg/mL)
Diethyl phthalate (0.25 mg/mL)	Di- <i>n</i> -octyl phthalate (0.5 mg/mL)
Dimethyl phthalate (0.25 mg/mL)	bis(2-Ethylhexyl)phthalate (0.5 mg/mL)

Method 8061A Phthalate Esters by GC/ECD

Phthalate Esters

M-8061-R1		1 x 1 mL
M-8061-R1-PAK	SAVE	5 x 1 mL
1.0 mg/mL each in Hexane		
bis(2- <i>n</i> -Butoxyethyl)phthalate	Dimethyl phthalate	
Butyl benzyl phthalate	Dinonyl phthalate	
Diamyl phthalate	Di- <i>n</i> -octyl phthalate	
Di- <i>n</i> -butyl phthalate	bis(2-Ethoxyethyl)phthalate	
Dicyclohexyl phthalate	bis(2-Ethylhexyl)phthalate	
Diethyl phthalate	bis(2-Methoxyethyl)phthalate	
Dihexyl phthalate	bis(4-Methyl-2-pentyl)phthalate	
Diisobutyl phthalate		

M-8061A		1 x 1 mL
M-8061A-PAK	SAVE	5 x 1 mL
1.0 mg/mL each in Hexane		

Butyl benzyl phthalate	Diethyl phthalate
bis(2-Ethylhexyl)phthalate	Dimethyl phthalate
Di- <i>n</i> -butyl phthalate	Di- <i>n</i> -octyl phthalate

Matrix Spike Solution

M-8061A-MS		1 x 1 mL
M-8061A-MS-PAK	SAVE	5 x 1 mL
2.0 mg/mL each in Acetone		

Butyl benzyl phthalate	bis(2-Ethylhexyl)phthalate
------------------------	----------------------------

Internal Standard

M-8061-IS		1 x 1 mL
M-8061-IS-PAK	SAVE	5 x 1 mL
5.0 mg/mL in Hexane		

Benzyl benzoate

Surrogate Standards

M-8061-SS		1 x 1 mL
M-8061-SS-PAK	SAVE	5 x 1 mL
50 µg/mL each in Acetone		
M-8061-SS-10X		1 x 1 mL
M-8061-SS-10X-PAK	SAVE	5 x 1 mL
500 µg/mL each in Acetone		

Dibenzyl phthalate	Diphenyl phthalate
Diphenyl isophthalate	



Thousands of Standards, just a click away

AccuStandard.com

Phthalates Solution Set

Set of 17 Phthalate Solutions
Each at 100 µg/mL Concentration

ALR-PHT-SET 17 x 1 mL
S in MeOH, S-CN in AcCN

Compound	CAS No.	Cat. No.	Unit
Benzyl butyl phthalate	85-68-7	ALR-082S	1 mL
Diamyl phthalate	131-18-0	ALR-098S	1 mL
Dicyclohexyl phthalate	84-61-7	ALR-099S	1 mL
Di(2-ethyl hexyl) phthalate (DEHP)	117-81-7	ALR-097S	1 mL
Diethyl phthalate	84-66-2	ALR-110S	1 mL
Di-hexyl phthalate	84-75-3	ALR-100S	1 mL
Diisodecyl phthalate	26761-40-0	ALR-101S	1 mL
Diisononyl phthalate	68515-48-0	ALR-102S	1 mL
Diisooctyl phthalate	27554-26-3	ALR-103S	1 mL
Dimethyl phthalate (DMP)	131-11-3	ALR-111S	1 mL
Di-n-butyl phthalate (DBP)	84-74-2	ALR-104S	1 mL
Di-n-octyl phthalate	117-84-0	ALR-105S	1 mL
Monobenzyl phthalate (mBzP)	2528-16-7	ALR-134S-CN	1 mL
Monobutyl phthalate (mBP)	131-70-4	ALR-135S-CN	1 mL
Monoethyl phthalate (mEP)	2306-33-4	ALR-137S-CN	1 mL
Monoethylhexyl phthalate (mEHP)	4376-20-9	ALR-138S-CN	1 mL
Monomethyl phthalate	4376-18-5	ALR-139S-CN	1 mL

Phthalate Esters Kit

PS-840C-R1-SET
Neats at 1 mL each.

15 vials

Calibration Mixture

PS-84C-1ML 1 mL vial
Neat at the stated weight %
Neat at stated weight %
Dimethyl phthalate 16.7%
Diethyl phthalate 16.7%
Di-n-propyl phthalate 33.3%
Di-n-butyl phthalate 33.3%

- (01) Dimethyl phthalate
- (02) Diethyl phthalate
- (03) Di-n-propyl phthalate
- (04) Di-iso-propyl phthalate
- (05) Di-n-butyl phthalate
- (06) Di-iso-butyl phthalate
- (07) Dipentyl phthalate
- (08) Dihexyl phthalate
- (09) Diheptyl phthalate
- (10) Dioctyl phthalate
- (11) Dinonyl phthalate
- (12) Didecyl phthalate
- (13) Diundecyl phthalate
- (14) Didodecyl phthalate

- (15) Phthalate Mixture PS-84C-1ML
Neat at stated weight %
Dimethyl phthalate 16.7%
Diethyl phthalate 16.7%
Di-n-propyl phthalate 33.3%
Di-n-butyl phthalate 33.3%

Miscellaneous Phthalate Mixes

Appendix IX Phthalate Mix

APP-9-PHTH-MIX

1000 µg/mL each in Cyclohexane

bis(2-Ethylhexyl)phthalate
Dibutyl phthalate
Di-n-octyl phthalate



Diisodecyl phthalate
Diisononyl phthalate
Benzyl butyl phthalate

1 x 1 mL
6 comps.

Phthalate Mix

ASM-146

1.0 mg/mL each in MeOH

Benzyl butyl phthalate
Dimethyl phthalate
Diethyl phthalate

Di-n-butyl phthalate
Di-n-octyl phthalate
bis(2-Ethylhexyl)phthalate

1 x 1 mL
6 comps.

Phthalate Esters Mix

M-PHE

M-PHE-PAK

At stated conc.(µg/mL) in Acetone

Benzyl butyl phthalate 10
bis(2-Ethylhexyl)phthalate 50
Diethyl phthalate 25

Dimethyl phthalate 25
Di-n-butyl phthalate 25
Di-n-octyl phthalate 50

1 x 1 mL
SAVE 5 x 1 mL
6 comps.



AccuStandard®

ISO Guide 34 ■ ISO/IEC 17025 ■ ISO 9001



CHROMATOGRAPHIC SPECIALTIES INC.

www.chromspec.com

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com

