



PFAS Workflow for *Solid Samples*

**X-TRACTION® PFAS (PFE*)
or MIX-TRACTION**
For Extraction



D-EVA Concentration
For the sensor controlled evaporation of extracts



FREESTYLE™ SPE PFAS
For clean-up up to 50 mL
(e.g. Dual SPE)



EluCLEAN® PFAS + **D-EVA Concentration**
SPE columns
For the sensor controlled evaporation to a few µL



* Pressurized fluid extraction



CHROMATOGRAPHIC **SPECIALTIES INC.**





PFAS Workflow for *Aqueous Samples*

FREESTYLE™ XANA PFAS

For clean-up of up to 4 L



FREESTYLE™ XANA PFAS TableTop


For clean-up of up to 250 mL



or

EluCLEAN® PFAS

SPE columns



D-EVA Concentration

For the sensor controlled evaporation to a few µL

