


### Essential Starter Kits with Push Caps and 2D Coded Tubes


A fully traceable storage solution at a low cost



SAVE  
up to  
30%!

### Starter Kits with Push Caps and Non-Coded Tubes

Economical solution for sample storage at temperatures down to -80°C



SAVE  
up to  
30%!

#### MJS BioLynx Order Information: MSMP50250C (0.75ml Tubes) / MSMP50251C (1.40ml Tubes)

0.75ml <b>or</b> 1.40ml 2D Data-Matrix Coded Screw Cap tubes V-bottom in Micronic 96-2 <b>or</b> 96-4 Rack with low cover	10 Racks
Pierceable TPE Capcluster Blue for capping 96 individual tubes	10 Mats
Handheld Wireless Scanner Mini	1 Unit
Traxcer Code Reader RS297	1 Unit
Manual Decapper-1	10 Pieces
Manual Decapper-8	10 Pieces

#### MJS BioLynx Order Information: MSMP50252C (0.75ml Tubes) / MSMP50253C (1.40ml Tubes)

0.75ml <b>or</b> 1.40ml Non-Coded Screw Cap tubes V-bottom in Micronic 96-2 <b>or</b> 96-4 Rack with low cover	10 Racks
Pierceable TPE Capcluster Blue for capping 96 individual tubes	10 Mats
Manual Decapper-1	10 Pieces
Manual Decapper-8	10 Pieces

### Essential Starter Kits with Screw Caps and 2D Coded Tubes

Convenient solution for maintaining a secure sample logistics system



### Starter Kits with Screw Caps and Non-Coded Tubes

Ideal for sample storage at ultra-low temperatures



MJS BioLynx Order Information: MSMP50254C (0.75ml Tubes) / MSMP50255C (1.40ml Tubes)	
0.75ml <b>or</b> 1.40ml 2D Data-Matrix Coded Screw Cap tubes V-bottom in Micronic 96-2 <b>or</b> 96-4 Rack with high cover	10 Racks
Screw Caps Blue	960 caps
Handheld Wireless Scanner Mini	1 Unit
Traxcer Code Reader RS297	1 Unit
Screw Cap Remover	10 Pieces

MJS BioLynx Order Information: MSMP50256C (0.75ml Tubes) / MSMP50257C (1.40ml Tubes)	
0.75ml <b>or</b> 1.40ml Non-Coded Screw Cap tubes V-bottom in Micronic 96-2 <b>or</b> 96-4 Rack with high cover	10 Racks
Screw Caps Blue	960 caps
Screw Cap Remover	10 Pieces