

MICREWTube®

Selection tool

With their unique models and incredible selection, Simport has developed the MicrewTube® Selection Tool to help you select the optimal product based on your needs.



Scan the QR code to access the MicrewTube® Selection Tool

Over 100 models to choose from.



MICREWTube® SELECTION

- Tubes are available from 0.5 mL to 2.0 mL.
- Offered sterile or non-sterile.
- Self-standing tubes or with conical bottoms.
- Sealing caps are available with washer seal or lip seal.
- Choose from a flat-top cap, with an attachment loop or without.
- CapInsert™ identification caps are also available.



MJS BioLynx Inc. is proud to represent the following companies in Canada ...

Cell & Molecular Biology

- Crystallography
- Immunohistochemistry
- Immunofluorescence
- Glycobiology
- Transfection & Electroporation
- 3D Cell Co-culture
- 3D Laboratory Services
- Density Gradient Media
- Cell Culture Media
- Cell Migration Assays
- Exosome Purification
- PCR, qPCR & PCR Clean Up Reagents



Antibodies and Proteins

- Enzymes & Antibodies
- Nucleosomes & Histones (Epigenetics)
- Lipids and Fatty Acids
- Recombinant Proteins
- ELISAs
- Protein Expression Reagents
- Cytokines



Labware & Equipment

- Pipette Tips
- Semi Automated Pipettors
- Freezer Racks
- Tubes
- Histology Labware
- Sample Collection and Management
- Sample Storage
- Microwell Plates
- Water Purification Systems
- Benchtop Equipment
- SPR
- Tube Labeling and Thawing
- qPCR



Microbiology

- Endotoxin Testing
- Glucan Products
- Yeast Growth Media
- Microscope Slides, Stains and Accessories
- Sterilization Products



+ More

- Detergents
- Protein Purification & Assay Reagents
- Chromatography



Please visit our website at www.biolynx.ca if you would like to be added to our mailing list.
1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca