

Purchase 1 Product and Get **50% OFF** a Second Product You've Never Tried Before!

Are there more Vector products you would like to try? Each time you purchase a Vector product between now and the end of November 2021, you can add on a Vector product you have never tried before and receive a 50% discount on that product.

Please contact us for assistance in choosing your new product.

Reference promotion code **VECT21NEW** at time of order.

Valid until November 30th, 2021.

50% discount is applicable to the list price and can only be applied to a product that has not been previously purchased by the customer.

Vector Laboratories offers a full range of products for:

- Immunohistochemistry / Immunocytochemistry
- Immunofluorescence
- Nucleic Acid Labeling and Detection
- Enzyme Substrates
- Lectins and Glycobiology Reagents
- Antibodies

Visit www.biolynx.ca/vector for more information.

Available in Canada from...

MJS
BioLynx
INC.

1-888-593-5969 • www.biolynx.ca • tech@biolynx.ca



FEATURED PRODUCTS



VECTASTAIN® Elite® ABC-HRP Kit (Peroxidase, Standard) **VECTPK6100**

The Elite® ABC series is approximately 5 times more sensitive than the original VECTASTAIN® ABC Kit with the same low background.

NEW! VectaMount® Express Mounting Medium

VECTH5700 (60mL)

Non-aqueous clearing and mounting medium enabling the rapid mounting of cell and tissue specimens following immunohistochemistry (IHC) staining. Saves time by eliminating the need for extensive ethanol and clearing washes prior to coverslipping.

NEW! VECTASHIELD® PLUS Antifade Mounting Media

VECTH190010 (10mL)

VECTH200010 - with DAPI (10mL)

VECTH19002 (2mL)

VECTH20002 - with DAPI (2mL)

New non-setting formulation for preserving fluorescence by inhibiting rapid photobleaching of fluorescent proteins and fluorescent dyes. Mounted sections can be viewed in 1 hour.

VECTASHIELD® Vibrance™ Antifade Mounting Media

VECTH17002 - no counterstain (10mL)

VECTH170010 - no counterstain (10mL)

VECTH18002 - with DAPI (2mL)

VECTH180010 - with DAPI (10mL)

Curing (setting) antifade mounting media for retention of Immunofluorescence staining on cell and tissue specimens. Visualize specimens in one (1) hour after mounting. No media retraction or blue background development on slides even with extended (several months) storage at room temperature.

ImmPACT® DAB EqV Substrate Kit, Peroxidase

VECTSK4103 (400mL)

VECTSK4103100 (100mL)

ImmPACT® DAB EqV is 3 - 5 times more sensitive than the original DAB substrate. Optimized for IHC/ICC/ISH. Heat Stable. Suitable for permanent mounting.

Vector® TrueVIEW™ Autofluorescence Quenching Kit

VECTSP8400 (15mL)

VECTSP8500 - with DAPI (15mL)

Quenching reagent that significantly reduces cell/tissue autofluorescence (AF) arising from Aldehyde fixation and/or tissue elements (collagen, elastin, RBCs).

M.O.M.® (Mouse on Mouse) Immunodetection Kit, Basic

VECTBMK2002

Significant reduction of endogenous mouse Ig staining when using mouse primary antibodies on mouse tissue. Useful in any mouse-on-mouse application, such as studies in genetically engineered mice, including transgenic and knock-out models.

H.O.H.™ (Human on Human) Immunodetection Kit

VECTHOH3000

The H.O.H. Immunodetection Kit is intended to detect human (or humanized) antibodies on frozen or paraffin embedded human tissue sections.

VECTASTAIN, Elite, VectaMount, VECTASHIELD, ImmPACT and M.O.M. are registered trademarks of Vector Laboratories, Inc.