


## 96 Round Well Royal Blue Sealing Mat with Spray Coated PTFE/Premium Silicone




229078



500246

This newest addition to our molded sealing mats is opaque for light sensitive solutions.

- Designed to fit our 96-Well Microplates with 7mm diameter round wells
- Silicone provides excellent resealability
- Resists coring and tearing after multiple needle injections
- Autoclavable
- Excellent chemical compatibility
- Manufactured from medical grade silicone
- The tight in-plate seal eliminates sample cross-contamination and evaporation
- Available pre-scored for easier needle penetration
- Alphanumeric grid references
- No coring or cross-talk
- Unscored mats are good with aggressive solvents and extreme shaking
- Pre-scored mats are good with aggressive solvents
- Available five (5) mats per case



Finneran sealing mats work in conjunction with the Porvair **Mat Capper** and **AutoCapper**. The Porvair **Mat Capper** is compact, portable, easy to use and requires only one operation to produce an accurate tight seal. The Porvair **AutoCapper** takes the effort out of applying friction sealing mats to deep well plates and tube racks. With the push of a button, the Porvair **AutoCapper** does the work, eliminating the manual application of mats to plates.

Part No.	Description
976057MR-96B	96 Round (7mm Diameter Plug) Royal Blue Sealing Mat with PTFE/Premium Silicone
976757MR-96B	96 Round Pre-Slit (7mm Diameter Plug) Royal Blue Sealing Mat with PTFE/Premium Silicone

### Mat Cappers for Friction Seals

Part No.	Description
229078	Mat Capper - Manual Cap Mat Applicator for storage and assay plates
500246	AutoCapper - Electrically operated Cap Mat Applicator for storage plates and vial racks, 110/220V

J.G. Finneran Associates, Inc. • Visit [www.JGFinneran.com](http://www.JGFinneran.com)  
Engineered and Manufactured in compliance with ISO 9001:2015 certified guidelines

Available in Canada from...

MJS  
**BioLynx**  
INC.

1-888-593-5969 • [www.biolynx.ca](http://www.biolynx.ca) • [tech@biolynx.ca](mailto:tech@biolynx.ca)

